

ANIMAL MARKER



The Animal Marker was developed for easy, long lasting identification of laboratory animals. It is intended for use on the fur, not on the skin.
(black animals excluded)

Main Features

- * Nontoxic and devoid of any hazardous materials
- * A single marker can be used to mark 300-500 rats (depending upon animal size and/or applied area)
- * Once applied, colors remain for 6-12 weeks.
- * Broad and Fine tip Markers available
- * Unlike permanent markers and art markers which can contain chemicals harmful to animals, Animal Marker use "legal pigments", which are approved and regulated by the Japanese Ministry of Health, Labor and Welfare, for use in food, drugs and cosmetics. These pigments are nontoxic and devoid of any hazardous materials.

CAUTION: Colors can stain skin and clothing. The use of gloves and lab coats is recommended.

New Item!

Animal Marker Black



Ordering Information

<Standard Tip>

Model №	Color
FG2000B	Blue
FG2000G	Green
FG2000R	Red
FG2000Y	Yellow
FG2000D	Orange
FG2000K	Black

NEW



Capacity	15 ml
Tip Size	15 x 10 x 10 mm
Minimum Order Quantity	2 ea (same color)/pack

<Fine Tip>

Model №	Color
FG2200B	Blue
FG2200G	Green
FG2200R	Red
FG2200Y	Yellow
FG2200D	Orange
FG2200K	Black


NEW



Capacity	10 ml
Tip Size	φ45 x 75 mm
Minimum Order Quantity	3 ea (same color)/pack

Specifications are subject to change without notice.

FG2200B
Blue used



MUROMACHI KIKAI CO., LTD.

4-2-1 Nihonbashi-Muromachi, Chuo-ku, Tokyo, 103-0022 JAPAN
E-mail:sales@muromachi.com

TEL:+81-3-3241-2444
FAX:+81-3-3241-2940