

LA PIAS LABORATORY ANIMAL EAR TAG



La Pias is a metal ear tag for small rodents and has following features :

(1) The Mouse ear tags come numbered from 1-100 (Available up to 999).
The Rat tags also come numbered from 1-100 (Available up to 999).

(2) 5 colors are available for both Mouse and Rat.



1. **Red** 2. **Green** 3. **Blue** 4. **Light Blue** 5. **Black**



Colored tags with laser etched numbers provide high visibility making them real easy to read.

(3) Once applied, they cannot be easily removed.

(4) Laser-etched number provides durability.

SPECIFICATIONS	
Material	Male Tag : Aluminum  Female Tag : Stainless Steel 
Size	Mouse : O.D.3.5mm with 3mm pin Rat : O.D.3.5mm with 5mm pin

PA-3 Ear Tag Applicator



Ordering Information

MODEL No.	DESCRIPTION	QUANTITY
M-3	La Pias Ear Tag for Mouse	100 ea/color (Numbered 1-100) or 100 ea/5colors (Numbered 1-20/color)
R-3	La Pias Ear Tag for Rat	100 ea/color (Numbered 1-100) or 100 ea/5colors (Numbered 1-20/color)
PA-3	Ear Tag Applicator for Mouse & Rat	1

NOTE : Minimum order quantity is 100 ea.

***Non colored Tag (Stainless Steel Tag) available.
Please feel free to contact us.**



Frequently Asked Questions

Question 1. Do La Pias ear tags fall out?

Answer: Although this is rare and typically only occurs when the ear tag back is not properly and/or tightly attached occasionally La Pias ear tags fall or are torn out. Ear tags are torn or fall out more frequently in group-housed mice in large part because of fighting and/or grooming behaviors.

Question 2. Can La Pias Ear Tags be used with neonatal pups?

Answer: La Pias Ear Tags are not generally suitable for use with neonatal pups as they tend to tear the ear.

Question 3. What is the difference between La Pias Aluminum Ear Tags and Stainless Steel Ear Tags?

Answer:

Material	color	Durability	Recommended Use
Aluminum	Five colors available (Red, Blue, Light Blue, Green, Black)	Only the Surface is hard making them less durable than stainless steel	Short-term-test-Single breeding
Stainless Steel	Silver only (Numbers in brown)	All parts are Hard and Durable	Short or long-term test Single or group breeding

Question 4. Are La Pias Ear Tags autoclavable?

Answer: Yes, La Pias Ear Tags are autoclavable; however, the ear tag applicator cannot be autoclaved.



Specifications are subject to change without notice.

MUROMACHI KIKAI CO., LTD.

4-2-1 Nihonbashi-Muromachi, Chuo-ku, Tokyo, 103-0022 JAPAN
E-mail:sales@muromachi.com

TEL:+81-3-3241-2444
FAX:+81-3-3241-2940

07 September, 2015