

MODEL MK-713/CU01 IGLOO FAST-TRACS RUNNING WHEEL with COUNTER



MK-713/CU01 IGLOO FAST-TRACS, Low-Cost Running Wheel is designed to measure spontaneous activity Levels of mice accurately. The plastic wheel is equipped with a sensor which detects the revolving of the wheel. As the wheel is mounted on the igloo (dome-shaped house) located in a home cage, it is easier for mice to run.

A fully automatic measurement can be done for multiple running wheels using a CompACT AMS Ver.3 Data Collection Software.

FTC-01 COUNTER UNIT



FTC-01 COUNTER UNIT
Display : LCD 6 digits (up to 999999)
Dimensions : W80 x D100 x H35 mm
Weight : 170 g

FT-713 RUNNING WHEEL	MC-713 MOUSE CAGE
Output : Contact	Water Bottle : 200 cc
Material : Plastic	Material : PC with Stainless Steel Lid
Dimensions : OD155 x H100 mm	Dimensions : W220 x D320 x H135 mm
Weight : Approx. 230 g	Weight : Approx. 2 kg

The Standard System Includes:	
FT-713 Igloo Fast Trak Running Wheel	1
MC-713 Mouse Home Cage	1
FTC-01 Counter Unit	1
Connection Cable	1
Water Bottle (200cc)	1

Specifications are subject to change without notice.

MUROMACHI KIKAI CO., LTD.

4-2-1 Nihonbashi-Muromachi, Chuo-ku, Tokyo, 103-0022 JAPAN
E-mail:sales@muromachi.com

TEL:+81-3-3241-2444
FAX:+81-3-3241-2940