

**Standard Pattern with Large Finger Loops**



80% of actual size



Straight  
Sharp/Blunt

14.5 cm  
Cutting edge: 27 mm

No. 14101-14

**Tungsten Carbide Mini-Lister**



For more information see page 211.



Angled to side  
Blunt/Blunt

9 cm  
Cutting edge: 15 mm

No. 14531-09



**Lister**

Excellent skin cutting scissors. Longer, blunt blade prevents damaging tissue underneath skin without losing the cutting line.



80% of actual size



Angled to side  
Blunt/Blunt

14 cm  
Cutting edge: 26 mm

No. 14131-14