


Moria Ultra Fine Forceps

These forceps have the most finely executed points available. They are extremely light, superbly balanced, and have a gentle spring action. The shanks are specially designed with an arch which ensures maximum force to the tips. The tips close flat for 6 mm from the points to give a large holding area. The tips have been expertly finished under magnification.



Moria MC40




Straight

11 cm

No. 11370-40

Moria MC40/B




Curved

11 cm

No. 11370-42



Moria 9980




Straight

13 cm

No. 11399-80

Moria 9987



Curved

13 cm


No. 11399-87

Moria Iris Forceps

Balanced with gentle spring action and tips that are perfectly matched.



Moria MC32/B


Smooth
Blunt

10 cm

Tip width: 0.5 mm

No. 11373-22

Moria MC31/B


Smooth
Blunt

10 cm

Tip width: 0.5 mm

No. 11373-12


Moria MC32



Microerrated

No. 11370-32

Moria MC31



Microerrated

No. 11370-31