


Insert Pins

All insert pins fit any of our pin holders.
Pins are antimagnetic and acid resistant.



Lancet Straight Curved
5 cm 5 cm 5 cm
Stainless steel Stainless steel Stainless steel
No. 26007-01 No. 26007-02 No. 26007-03



Lancet half Spatula
5 cm 5 cm
Stainless steel Stainless steel
No. 26007-04 No. 26007-05

Black Anodized Minutien Pins

These high quality pins are black for best visibility
under high magnification and bright light.

Tip diameter: 0.0175 mm
Rod diameter: 0.15 mm
12 mm

Box of 500 pins
No. 26002-15

**MOST
POPULAR**



Stainless Steel Minutien Pins

The thinnest and finest pins available.

**NEW
ITEM**

Tip diameter: 0.0125 mm
Rod diameter: 0.1 mm
12 mm

Box of 500 pins
No. 26002-10

**MOST
POPULAR**

Tip diameter: 0.0175 mm
Rod diameter: 0.15 mm
12 mm

Box of 500 pins
No. 26003-15

**MOST
POPULAR**


Tip diameter: 0.02 mm
Rod diameter: 0.2 mm
12 mm

Box of 500 pins
No. 26002-20

**MOST
POPULAR**

Black Anodized Insect Pins

Nylon heads
Tip diameter: 0.03 mm
Pkg. of 100 pins



Size	000	00	0	1	2	3	4	5	6	7
Rod diameter (mm)	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.65	0.70
No. 26000	-25	-30	-35	-40	-45	-50	-55	-60	-65	-70


Stainless Steel Insect Pins

Nylon heads
Tip diameter: 0.03 mm
Pkg. of 100 pins

Size	000	00	0	1	2	3	4	5	6	7
Rod diameter (mm)	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.65	0.70
No. 26001	-25	-30	-35	-40	-45	-50	-55	-60	-65	-70

Etched Tungsten Micro Needles

These tungsten micro needles are drawn to a 1 micron tip and three different rod diameters are offered. Each provides lateral stability as well as minimal spring action during manipulation. The micro needles can be used with our pin holders shown on page 114.



Tip diameter: 1 μ m (0.001 mm)
12 mm

Box of 10 needles

Rod diameter (mm)	Item	Price
0.125	No. 10130-05	
0.25	No. 10130-10	
0.5	No. 10130-20	