

Alm Retractor with Wire Teeth



7 cm
Maximum spread: 4 cm
Arms gently curved down


No. 17010-07

Alm Retractor with Blunt Teeth



7 cm
Maximum spread: 4 cm
Arms gently curved down

Blunt teeth
No. 17008-07



7 cm
Maximum spread: 4 cm
Arms gently curved down

Long blunt teeth
No. 17009-07

Alm Retractor with Blades



7.5 cm
Maximum spread: 4.5 cm
Arms gently curved down
Thin, curved blades

No. 17009-08