

Vannas-Tübingen



15003-08



**Straight  
Sharp/Sharp**

8.5 cm  
Cutting edge: 5 mm  
Tip Ø: 0.075 mm

**No. 15003-08**



**Angled up  
Sharp/Sharp**

9.5 cm  
Cutting edge: 5 mm  
Tip Ø: 0.075 mm

**No. 15005-08**



**Curved  
Sharp/Sharp**

8.5 cm  
Cutting edge: 5 mm  
Tip Ø: 0.075 mm

**No. 15004-08**



**Angled to side  
Sharp/Sharp**

8.5 cm  
Cutting edge: 5 mm  
Tip Ø: 0.075 mm

**No. 15008-08**



**Straight  
Sharp/Sharp**

8.5 cm  
Cutting edge: 5 mm  
Tip Ø: 0.075 mm

**No. 15903-08**



**Curved  
Sharp/Sharp**

8.5 cm  
Cutting edge: 5 mm  
Tip Ø: 0.075 mm

**No. 15904-08**

Vannas Spring Scissors High Polished

The highly polished surface makes the scissors more resistant to corrosion and easier to clean. It also prevents tissue from sticking on it while operating, it's also called non-sticky surface.



Vannas Spring Scissors Round Handle

The round handle allows an easy finger adjustment during the application. The hollow grind in the handle reduces the weight. The fine blades are tapered to give a better overview.



**Straight  
Sharp/Sharp**

9 cm  
Cutting edge: 5 mm  
Tip Ø: 0.075 mm

**No. 15403-08**



**Curved  
Sharp/Sharp**

9 cm  
Cutting edge: 5 mm  
Tip Ø: 0.075 mm

**No. 15404-08**