

**Noyes**



**Straight  
Sharp/Sharp**

12 cm  
Cutting edge: 14 mm  
Tip Ø: 0.2 mm


**No. 15012-12**



**Curved  
Sharp/Sharp**

12 cm  
Cutting edge: 14 mm  
Tip Ø: 0.2 mm

**No. 15011-12**



**Angled up  
Sharp/Sharp**

12 cm  
Cutting edge: 14 mm  
Tip Ø: 0.2 mm

**No. 15013-12**

**Tungsten Carbide Noyes**

Tungsten carbide cutting edges provide long lasting sharpness.



**Straight  
Sharp/Sharp**

12 cm  
Cutting edge: 14 mm  
Tip Ø: 0.275 mm

**No. 15514-12**

**Castroviejo**



**Curved  
Sharp/Sharp**

10 cm  
Cutting edge: 10 mm  
Tip Ø: 0.2 mm

**No. 15017-10**