

**Standard Pattern with Large Finger Loops**



80% of actual size



**Straight Sharp/Blunt**

14.5 cm  
Cutting edge: 30 mm

**No. 14101-14**

**Tungsten Carbide Mini-Lister**



**Angled to side Blunt/Blunt**

9 cm  
Cutting edge: 14 mm

**No. 14531-09**

**Lister**

Excellent skin cutting scissors. Longer, blunt blade prevents damaging tissue underneath skin without losing the cutting line.



80% of actual size



**Angled to side Blunt/Blunt**

14 cm  
Cutting edge: 26 mm

**No. 14131-14**