


Friedmann-Pearson


This popular, fine-tipped rongeur is suited for delicate work. The fine-cupped jaws have been designed with a longer and slimmer approach for better visibility.



Cautionary note: The fine rongeur tips shown on this page are extremely fragile. Removal of material with a twisting motion can damage the tips.



200%



Straight

14 cm
Cup size: 0.5 mm

No. 16220-14

Curved

14 cm
Cup size: 0.5 mm

No. 16221-14

14 cm
Cup size: 0.7 mm

No. 16120-14

14 cm
Cup size: 0.7 mm

No. 16121-14

14 cm
Cup size: 1 mm

No. 16020-14

14 cm
Cup size: 1 mm

No. 16021-14

