

**Bone Forceps**

Strong cross-serrated tips for bone crushing and removal.



**Angled**

11.5 cm  
Tip Ø: 1 mm  
Cross serrated

**No. 16060-11**

**Miniature Self-Tapping Bone Screws**

Used extensively to anchor dental cement to bone, these self-tapping screws have superior holding strength over ordinary machine screws. Stainless steel with fillister flat head.



4 mm  
Shaft diameter: 0.86 mm  
Pkg. of 100 screws

**No. 19010-10**

4.8 mm  
Shaft diameter: 1.19 mm  
Pkg. of 100 screws

**No. 19010-00**

Pkg. of 10 screws

**No. 19010-11**

Pkg. of 10 screws

**No. 19010-01**

For pre-drilling use our 0.7 mm burrs, see page 181.

For pre-drilling use our 0.9 mm burrs, see page 181.

**Bone Wax**

A beeswax-based compound for bone hemostasis. Soft, kneadable and ready to use. **Sterile.**



Pkg. of 12 foil packs (2,5 g each)

**No. 19009-00**