

SCALPELS

**Carbon Steel Scalpel**

Cold forged from carbon steel and sharpened to eliminate tissue distortion and blade drag. 12.5 cm plastic handle with slip-proof grip. **Sterile.**



10316-14  
10315-12



**22.5° Cutting angle**

13.5 cm  
Cutting edge: 4.5 mm  
Tip thickness: 0.15 mm


**No. 10316-14**



**15° Cutting Angle**

13 cm  
Cutting edge: 5 mm  
Tip thickness: 0.15 mm

**No. 10315-12**



**Crescent blade**

14 cm  
Cutting edge: 11 mm  
Tip thickness: 0.2 mm

**No. 10317-14**



**Needle blade**

14.5 cm  
Cutting edge: 4 mm  
Tip thickness: 0.2 mm

**No. 10318-14**

**Interchangeable Carbon Steel Blades**



Handle for Interchangeable Blades  
9.5 cm

**No. 10035-00**

**No. 10035-12**



**Angled Carbon steel**

Cutting edge: 10 mm  
Pkg. of 5 blades

**No. 10035-15**



**Curved Carbon steel**

Cutting edge: 10 mm  
Pkg. of 5 blades

**No. 10035-20**



**Small / Flat Carbon steel**

Cutting edge: 1.5 mm  
Pkg. of 5 blades

**No. 10035-05**



**Large / flat Carbon steel**

Cutting edge: 2.75 mm  
Pkg. of 5 blades

**No. 10035-10**



**Cane Carbon steel**

Cutting edge: 12 mm  
Pkg. of 5 blades

**No. 10035-25**



**Arc Carbon steel**

Cutting edge: 12 mm  
Pkg. of 5 blades

**No. 10035-30**