

**Scalpel Handle #4 with Silicone**

For heavy, gross anatomy work.  
Molded to fit hand.



13.5 cm

No. 10009-13

**Standard Scalpel Handle #4 Solid**

Especially suitable for large dissections.



13 cm

No. 10004-13

**Blades for #4 Scalpel Handle**

The alloy used has a high carbon content which provides a longer lasting edge. Blades come individually wrapped. **Non-sterile.**



Box. of 100 blades

No. 10020-00



Box. of 100 blades

No. 10021-00



Box. of 100 blades

No. 10022-00



Box. of 100 blades

No. 10023-00