


**Dumont #55**

Bevelled, like #5 but with extra light and fine shanks.



**Biology**

11 cm  
Inox  
0.05 x 0.02 mm


**No. 11255-20**

11 cm  
Dumostar  
0.05 x 0.02 mm

**No. 11295-51**

**Dumont #5  
Mirror Finish Forceps**


The last 5 mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves or any other adhesive tissue.



**Standard**

11 cm  
Inox  
0.1 x 0.06 mm

**No. 11251-23**



**Biology**

11 cm  
Inox  
0.05 x 0.02 mm

**No. 11252-23**

**Dumont #5**

These are the most widely used Dumont forceps.  
Bevelled, very fine shanks for best visibility.




11251-20  
11252-20



11252-40



11251-10



**Standard**

11 cm  
Inox  
0.1 x 0.06 mm


**No. 11251-20**

Dumoxel

**No. 11251-30**

Dumostar

**No. 11295-00**



**Standard**

11 cm  
Carbon steel  
0.08 x 0.04 mm

**No. 11251-10**



**Biology**

11 cm  
Inox  
0.05 x 0.02 mm

**No. 11252-20**

Dumoxel

**No. 11252-30**

Dumostar

**No. 11295-10**

Titanium

**No. 11252-40**