

Dumont WA #5


These specially designed forceps offer improved stability with an ergonomic body that allows the forceps to fit more comfortably in your hand. A domed end with no soldering grooves allows the forceps to be thoroughly and easily cleaned, preventing biofilm residues.

Dumont #5XL

FORCEPS/TWEEZERS



Side View



Standard

12 cm
Inox
0.10 x 0.06 mm


No. 11251-21



Biology

12 cm
Inox
0.05 x 0.02 mm

No. 11252-21



Standard

15 cm
Inox
0.1 x 0.06 mm

No. 11253-10