


Dumont #L5

Slide sleeve keeps tips together so specimens may be securely held and retained. Slide may be locked in place so that forceps can be used in regular fashion.




Standard

11 cm
Inox
0.1 x 0.06 mm

No. 11253-27

**Dumont #5
Teflon® Coated**

These antimagnetic forceps are ideal for use in cryogenic environments and for EM samples. Teflon® coating provides a non-stick surface.




Standard

11 cm
Teflon coated Inox
0.1 x 0.06 mm

No. 11626-11

Dumont #5/45

45° tips for improved visibility and control.



Standard

11 cm
Dumoxel
0.1 x 0.06 mm

No. 11251-35

**Dumont #5/45
Cover Slip Forceps**

Designed to pick up cover slips out of multi-well plates. The tip width is squared for a better grip.



300%



1 mm from tip to bend

11.5 cm
Dumoxel
Upper tip: 0.1 x 0.08 mm
Lower tip: 0.1 x 0.02 mm

No. 11251-33

FORCEPS/TWEEZERS