


Dumont #4

Similar to #5, but with broader shanks.




Standard

11 cm
Dumoxel
0.13 x 0.08 mm

No. 11241-30

Dumostar

No. 11294-00




Biology

11 cm
Dumoxel
0.06 x 0.04 mm

No. 11242-40

**Dumont #2
Laminectomy Forceps**

These forceps can be used to break off delicate bone during mouse laminectomies, as well as fine skull bone removal.




Standard

12 cm
Inox
0.4 x 0.02 mm

No. 11223-20

**Dumont #2AP
Epoxy Coated**

Provides good insulation for up to 220 AC voltage work and highly impervious to chemical solutions. Frozen or super cold sections will not adhere to the gripping surfaces.



Standard

12 cm
Epoxy coated Inox
1.5 x 0.7 mm

No. 11220-21