

Gold Plated Forceps

Gold-plating gives corrosion resistance to most chemicals, salts and acids. These forceps are perfect for any application where corrosion resistance and toughness are primary requirements.



Straight

11 cm
Gold plated
Stainless Steel
0.1 x 0.06 mm

No. 11915-11



Curved

12 cm
Gold plated
Stainless Steel
0.2 x 0.1 mm

No. 11917-12




Angled 30°

11 cm
Gold plated
Stainless Steel
0.1 x 0.06 mm

No. 11911-11

Diamond Coated Forceps


Diamond coating provides increased corrosion resistance and greater protection from wear and tear, which extends the operational life of the forceps, as well as providing a stronger grip at the tips.



Straight

14 cm
Diamond coated
Stainless Steel
0.15 x 0.1 mm

No. 11990-14



Angled 30°

14 cm
Diamond coated
Stainless Steel
0.15 x 0.1 mm

No. 11991-14



Straight

12 cm
Diamond coated
Stainless Steel
1.65 x 0.05 mm

No. 11992-12

FORCEPS/TWEEZERS