


## FINE FORCEPS

### Mirror Finish Broad Shanks

These forceps have long and narrow shanks with thin, highly polished tips to prevent tissue from sticking. Ideal for embryos, nerves, or other adhesive tissues.



**Straight**

13 cm  
0.3 x 0.15 mm

No. 11923-13


### #2A Forceps



**Straight**

12 cm  
1.5 x 0.13 mm

No. 11992-20




**Curved**

13 cm  
1.5 x 0.1 mm

No. 11972-20

### Forceps with Replaceable Tips

Stainless steel body with replaceable plastic tips. Tips are made of polyvinylidene fluoride (PVDF) and are electrostatic discharge (ESD) safe. Body and plastic tips are corrosion and abrasion resistant.



**Straight**

13 cm  
0.5 x 0.5 mm

Stainless Steel body  
with plastic tips

No. 11980-13



Pkg. of 2 tips and 3 screws  
Replacement plastic tips

Straight / Fine

4 cm  
0.6 x 0.6 mm

No. 11988-01



Pkg. of 2 tips and 3 screws  
Replacement plastic tips

Straight / Very fine

4 cm  
0.5 x 0.5 mm

No. 11988-02



Pkg. of 2 tips and 3 screws  
Replacement plastic tips

Curved / Very fine

4 cm  
0.5 x 0.5 mm

No. 11988-03