


Fine Forceps Broad Shanks

Straight forceps with sturdy coarse tips.




Straight

12 cm
0.2 x 0.12 mm

No. 11445-12

Fine Forceps Mirror Finish


The last 5 mm on the inside of the tips are highly polished to help prevent tissue from sticking. Ideal for embryos, nerves or other adhesive tissue.



Straight

11 cm
0.1 x 0.06 mm

No. 11412-11




Straight

11 cm
0.05 x 0.02 mm

No. 11413-11

Fine Forceps Self-Closing


Ideal for lifting and holding EM grids. Consistent pressure at tips. Correctly tensioned "cross-over" shanks come to a microscopically matched point



Straight

11.5 cm
0.08 x 0.04 mm

No. 11480-11



Straight

11 cm
0.1 x 0.06 mm

No. 11485-11



Curved

11.5 cm
0.2 x 0.12 mm

No. 11487-11

FORCEPS/TWEEZERS