



**Bonn**



20x magnification




**Straight**

Smooth with platform

7 cm

Tip width: 0.3 mm  
Tying platform: 5 mm

No. 11083-07



**Straight**

1x2 Teeth


7 cm

Tip width: 0.3 mm  
Teeth width: 0.1 mm

No. 11084-07

**Mirror Finish Bonn Forceps**

The fine Bonn forceps features a tying platform and is also available with additional 1x2 teeth. It's designed for tissue manipulation and suturing. The wide handle ensures safe handling. The fine and precise tips make it suitable for work under a high magnification.



20x magnification




**Straight**

Tying platform with 1x2 teeth

10.5 cm

Tip width: 0.12 mm  
Tying platform: 5 mm

No. 11090-10



**Straight**

Smooth with platform


10.5 cm

Tip width: 0.3 mm  
Tying platform: 5 mm

No. 11093-10

**Mirror Finish Suture Tying Forceps with Teeth**

A unique tip style fusing a suture tying platform and a 1x2 teeth tip for tissue manipulation.



20x magnification



**Straight**

Tying platform with 1x2 teeth

9 cm

Tip width: 0.3 mm  
Tying platform: 6 mm

No. 11090-09

FORCEPS/TWEEZERS