

Moria Suture Tying Forceps

Used to suture after surgery is complete. This requires the manipulation of a calibrated nylon thread and the use of a specific instrument: a tying forceps, straight or curved, without teeth.



Moria 9604




Straight

Smooth with platform

9.5 cm
Tip Ø: 0.3 mm
Tying platform: 8 mm

No. 11393-04

Moria 9605



Curved

Smooth with platform

8.5 cm
Tip Ø: 0.3 mm
Tying platform: 7.5 mm

No. 11393-05

McPherson Forceps



Angled 45°


Smooth with platform

9 cm
Tip Ø: 0.4 mm
Tying platform: 6 mm

No. 11063-07

Suture Tying Forceps

The platform tips close flat. This ensures a reliable grip on sutures.



Straight

Smooth with platform

13.5 cm
Tip Ø: 0.3 mm
Tying platform: 6 mm

No. 00272-13

Suture Tying Forceps

The platform tips close flat. This ensures a reliable grip on sutures.



Straight

Smooth with platform

10 cm
Tip Ø: 0.5 mm
Tying platform: 5 mm

No. 18025-10