

Plastic Forceps

These forceps are non-magnetic, acid resistant and able to withstand temperatures up to 240°C (464°F). Tips do not spread under pressure. Plastic forceps are made of Acetyl and will be sold in packages of 10 pieces.



Straight
Sharp

11.5 cm
0.4 x 0.25 mm

No. 11700-00



Straight
Blunt

11.5 cm
1.5 x 0.4 mm

No. 11700-03



Curved
Sharp

11.5 cm
0.3 x 0.2 mm

No. 11700-01



Angled 45°
Sharp

11.5 cm
0.3 x 0.2 mm

No. 11700-02




Curved
Blunt

11.5 cm
1.5 x 0.4 mm

No. 11700-04

Cross Action Forceps

These tissue transport forceps have a gentle action.



Angled
Cross serrated

10 cm
4 x 8 mm

No. 11075-00

Cover Glass Forceps



Straight
Smooth

10.5 cm
4 x 10 x 0.6 mm

No. 11074-02



Angled
Smooth

10.5 cm
4 x 10 x 0.6 mm

No. 11073-10

FORCEPS/TWEEZERS