

**Ring Forceps**

For transport of tissue and organs. With limit screw to prevent excessive applied pressure.

**Circle Forceps**

Cupped forceps for atraumatic gripping of organs and tissues.

**Live Insect Handling Forceps**

Very delicate.

FORCEPS/TWEEZERS



**Straight**  
Ring

8.5 cm  
ID: 1 mm; OD: 1.5 mm  
x 0.5 mm

No. 11101-09



**Straight**  
Ring

9 cm  
ID: 2.2 mm; OD: 3 mm  
x 1.3 mm

No. 11103-09



**Straight**  
Ring

9 cm  
ID: 4.8 mm; OD: 6 mm  
x 1.5 mm

No. 11106-09



**Straight**  
Circle

11.5 cm  
7 x 0.5 mm

No. 11072-11



**Straight**  
Narrow

10 cm  
0.9 x 0.15 mm

No. 26029-10



**Straight**  
Wide

10 cm  
2.3 x 0.3 mm

No. 26030-10