


Adson DeBakey

Atraumatic DeBakey Serrations

This tooth pattern, as shown, consists of one row of atraumatic teeth meshed with two opposite rows of atraumatic teeth.

These two forceps are ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures.



Straight
DeBakey serrations

12.5 cm
1.75 x 1.5 mm

No. 11617-12



Straight
DeBakey serrations

15 cm
1.75 x 1.5 mm

No. 11616-15

Adson-Graefe




Straight
Fine teeth

12.5 cm
4.5 x 1.25 mm

No. 11030-12

Adson-Brown




Straight
Shark teeth

12 cm
2 x 1.4 mm

No. 11627-12