

Semken

The Semken tweezers are characterized by their slim shape and fine tips. Their anatomical (toothing) or surgical (1x2 teeth) tips make them ideal for grasping tissue. The serrated grip surface ensures a secure hold during use.



11608-15



FORCEPS/TWEEZERS



Straight
Serrated

13 cm
Tip width: 1.3 mm

No. 11008-13



Curved
Serrated

13 cm
Tip width: 1.3 mm

No. 11009-13



Straight
1x2 Teeth

12.5 mm
Tip width: 1.3 mm

No. 11024-12

15 cm
Tip width: 1.5 mm
Titanium

No. 11608-15

15 cm
Tip width: 1.5 mm

No. 11009-15

DeBakey Round Handle Forceps

Ideal for gentle, atraumatic gripping of vessels, tissue and other delicate structures. Two rows of atraumatic teeth interfit with one opposite row of atraumatic teeth. The round handles allow easy fingertip adjustments.



Straight
DeBakey serrations

15 cm
Tip width: 1.2 mm

No. 11076-15