

## CLAMPS

### Micro Clip Applicator with Lock



14 cm

For use with FST 18055 and 18555 series micro serrefines.

Lock: Yes

**No. 18056-14**

### Forceps Style Clip Applicator

This applicator features special U-shaped side supports for improved seating of the clips.



14 cm

For use with FST 18055 and 18555 series micro serrefines.

Lock: No

**No. 18057-14**

### Biemer Clip Applicator (Forceps Style)

Application forceps for Biemer clamps made of surgical stainless steel.



14 cm

For use with FST 18058 (Biemer Clip) series

Lock: No

**No. 18059-14**

### Stainless Steel Micro Serrefines

Used to obtain a clear surgical field for extended arterial cross-clamping.



Micro Serrefines should never be compressed with fingers; use only the appropriate applicator No. 18056-14 or No. 18057-14.

Delicate atraumatic serration



Jaw dimensions (mm)	Total length (mm)	Jaw pressure (g)	Item
10 x 2	19	60	<b>No. 18055-01</b>
8 x 2	16	85	<b>No. 18055-02</b>
6 x 1	15	100	<b>No. 18055-03</b>
4 x 0.75	13	125	<b>No. 18055-04</b>

Jaw dimensions (mm)	Total length (mm)	Jaw pressure (g)	Item
6 x 1	17	85	<b>No. 18055-05</b>
4 x 0.75	16	85	<b>No. 18055-06</b>

### Titanium Micro Serrefines

Lightweight with strong clamping pressure.

Delicate atraumatic serration



Jaw dimensions (mm)	Total length (mm)	Jaw pressure (g)	Item
10 x 2	19	60	<b>No. 18555-01</b>
8 x 2	16	85	<b>No. 18555-02</b>
6 x 1	15	100	<b>No. 18555-03</b>
4 x 0.75	13	125	<b>No. 18555-04</b>