



**Moria Vessel Clamp**

Straight  
33 mm  
Serrated: Yes  
Light pressure  
Jaw dimensions: 13 x 1.6 mm



Moria 201/A  
**No. 18320-11**

**Moria Bulldog Type Clamps**

53 mm  
Serrated: No  
Strong pressure  
Jaw dimensions: 16 x 2.5 mm



Straight  
Moria MC43  
**No. 18374-43**



Curved  
Moria MC44  
**No. 18374-44**

**Bulldog Artery Spring Clamp**

Can be used with a minimum of distortion for MRI investigations.



Straight  
30 mm  
Serrated: Yes  
Light pressure



**No. 18144-30**

**Dietrich Bulldog Clamps**

Serrated: Yes

Closing pressure: 50 g



Straight  
52.5 mm  
**No. 18038-45**



Angled  
50 mm  
**No. 18039-45**

**Bulldog Type Serrefines**

Serrated: Yes

Straight  
28 mm  
Medium pressure



**No. 18050-28**

Curved  
28 mm  
Medium pressure



**No. 18051-28**

Straight  
35 mm  
Strong pressure



**No. 18050-35**

Straight  
50 mm  
Strong pressure



Curved  
35 mm  
Strong pressure



**No. 18051-35**

Curved  
50 mm  
Strong pressure

**No. 18051-50**

**Schwartz Micro Serrefines**

Serrated: No  
Strong pressure  
26 mm



Straight  
**No. 18052-01**



Light bend  
**No. 18052-02**



Sharp bend  
**No. 18052-03**



**Biemer Clip**

Serrated: Yes  
Straight

Closing pressure: 0.29 - 0.39 N / 30 - 40 g  
15 mm

Clamping Length: 6 mm

**No. 18058-15**



18 mm  
Clamping Length: 9 mm

**No. 18058-18**



**Mini Serrefines**

28 mm  
Medium pressure  
Serrated: Yes  
Jaw dimensions: 6 x 2 mm

Straight  
**No. 18054-28**



Curved  
**No. 18053-28**



VASCULAR INSTRUMENTS