

Vessel Dilators

Highly modified forceps with elongated, parallel tips which are highly polished. Designed for controlled intraluminal vessel dilation.



**Angled 10°**

11 cm  
Tip Ø: 0.1 mm

**No. 00126-11**

Tip Ø: 0.2 mm


**No. 00125-11**



**Angled 45°**

13.5 cm  
Tip Ø: 0.2 mm

**No. 00276-13**



**Angled 10°**

12 cm  
Tip Ø: 0.1 mm

**No. 18131-12**

Tip Ø: 0.2 mm

**No. 18132-12**