


**V-Hook**

Aids in feeding fine suture around blood vessels, works very similar to ligation aids.



300%



11 cm  
Tip Dimensions: 0,4 x 0,2 mm

**No. 18052-12**

**Ligation Aids**

A handy instrument for leading a fine suture under and around a blood vessel for tying-off purposes. Also useful for ligating peripheral nerves. Grooves toward hole help threading.



300%



10.5 cm  
Hole Ø: 0.4 mm

**No. 18062-12**



200%

200%



13.5 cm  
Hole Ø: 1 mm

**No. 18062-13**



13.5 cm  
Hole Ø: 1 mm

**No. 18063-13**