


Vessel Cannulation Forceps

These vessel cannulation forceps are useful whenever a fine plastic tube has to be introduced into a small blood vessel of almost equal size. The hollowed jaws hold the tubing securely without deforming it. The tip of the tubing is directed exactly into the vessel opening without damage to the vessel from unwanted movement.



11 cm
For tubing of 0.35 mm OD
No. 00608-11

11 cm
For tubing of 0.5-1 mm OD
No. 00574-11



13 cm
For tubing of 0.35 mm OD
No. 11280-11

11 cm
For tubing of 0.5 mm OD
No. 11281-11

11 cm
For tubing of 1 mm OD
No. 11282-11



11.5 cm
For tubing of 0.35 mm OD
No. 18403-11

11.5 cm
For tubing of 0.6 mm OD
No. 18406-11

13.5 cm
For tubing of 0.6 mm OD
No. 18406-13