

**Castroviejo**

Micro pattern, with slightly blunted tips.



**Straight**

9 cm  
Lock: Yes

**No. 12060-01**

Lock: No

**No. 12060-02**



**Curved**

9 cm  
Lock: Yes

**No. 12061-01**

Lock: No

**No. 12061-02**

**Kalt Needle Holder**

The Kalt needle holder is characterized by its extraordinary design and good handling. The closure system closes softly and securely and can be easily opened again with your thumb.



**Straight**

Serrated: Yes  
13.5 cm  
Lock: Yes

**No. 12080-14**

**Curved**

Serrated: Yes  
13.5 cm  
Lock: Yes

**No. 12081-14**



Never use a micro needle holder with needles that are too thick, as this can easily damage the fine jaws. Use only micro needles of less than 0.2 mm diameter (see page 148).