


**Suture Needles**

- Polished stainless steel
  - 3/8 curved
  - Round with standard eye
  - Suitable for sutures up to 0.32 mm diameter
- Ø: 0.35 mm

**Sterile: No**

Pkg. of 12 needles



Size	Item
00	No. 12050-20
0	No. 12050-10
1	No. 12050-01
2	No. 12050-02
3	No. 12050-03

**Silk Suture**

Silk suture has a smooth flow through tissue while maintaining knot security.

Braided Silk

**Sterile: No**

Can be ETO sterilized



Size	Suture Ø (mm)	Spool length (m)	Item
7/0	0.05	91	No. 18020-70
6/0	0.08	23	No. 18020-60
5/0	0.12	23	No. 18020-50
4/0	0.17	23	No. 18020-40
3/0	0.22	23	No. 18020-30
2/0	0.32	23	No. 18020-20
0	0.37	23	No. 18020-00
1	0.45	23	No. 18020-01
2	0.55	23	No. 18020-02
3	0.65	23	No. 18020-03

WOUND CLOSURE

For Suture Tying Forceps, see pages 79, 80 and 81.

**Nylon Monofilament Suture**

Nylon is known for its uniform diameter flexibility and outstanding consistent knot security.

Nylon monofilament

**Non sterile**

Can be ETO sterilized



Size	Suture Ø (mm)	Spool length (m)	Item
8/0	0.04	91	No. 18030-80
7/0	0.05	91	No. 18030-70
6/0	0.07	91	No. 18030-60
5/0	0.1	91	No. 18030-50
4/0	0.15	91	No. 18030-40
3/0	0.2	91	No. 18030-30