

Reusable Gavage Feeding Needles

Stainless steel feeding needles with luer lock hub. Also ideal for transcatheter perfusions with a fixative. The needles come with a round tip.



Round tip

No. 18061-24



Round tip

No. 18061-50

Gauge	Length (mm)	Tip Ø (mm)	Tip Shape	Shape	Mice (g)	Rats (g)	Guinea pigs (g)	Item
24	25	1.25	Round	Straight	>14			No. 18061-24
22	25	1.25	Round	Straight	15-20			No. 18061-22
20	38	2.25	Round	Straight	20-25	50-75		No. 18064-20
20	38	2.25	Round	Curved	20-25	50-75		No. 18065-20
18	38	2.25	Round	Straight	20-30	100-200	250-300	No. 18061-38
18	50	2.25	Round	Curved	25-30	100-200	250-300	No. 18061-50
16	75	3	Round	Curved		150-300	350-450	No. 18061-75
16	100	3	Round	Curved		150-300	350-450	No. 18061-10

These needles come with a barrel tip.



Barrel tip

No. 18060-20



Gauge	Length (mm)	Tip Shape	Shape	Mice (g)	Rats (g)	Guinea pigs (g)	Item
24	25	Barrel	Straight	>14			No. 18060-24
22	25	Barrel	Straight	15-20			No. 18060-22
20	30	Barrel	Straight	20-25	50-75		No. 18060-20
20	38	Barrel	Straight	20-25	50-75		No. 18063-20
20	38	Barrel	Curved	20-25	50-75		No. 18067-20
18	38	Barrel	Straight	20-30	100-200	250-300	No. 18060-18
18	50	Barrel	Curved	25-30	100-200	250-300	No. 18067-50
16	75	Barrel	Curved		150-300	350-450	No. 18067-75
16	100	Barrel	Curved		150-300	350-450	No. 18067-10