

## EAR PUNCHES

### Ear Punch with Replaceable Tip

Designed to ensure a fast and clean punch, this ear punch is lightweight and fits comfortably in the palm of your hand. The tip is replaceable. **Autoclavable.**



#### Handle and Tip

11.5 cm  
Anodized aluminum handle  
and stainless steel tip  
Punch Ø: 2 mm

**No. 24213-02**

#### Replacement Tip

Punch Ø: 2 mm

**No. 24213-12**

### Ear Punch

Designed to ensure a fast and clean punch, this ear punch is lightweight and fits comfortably in the palm of your hand. FST Exclusive Item. **Autoclavable.**



11.5 cm  
Stainless Steel and anodized  
aluminum  
Punch Ø: 2 mm

**No. 24214-02**

### Finger Loop Ear Punch

This ear punch provides stability and control for researchers using the numbering system. The 1 mm ear punch offers greater numbering variations and/or patterns.



10.5 cm  
Punch Ø: 2 mm

**No. 24212-02**

Punch Ø: 1 mm

**No. 24212-01**