

## INSTRUMENT STORAGE

### Roll Up Travel Case

Vinyl dissection roll up case with Velcro strap to hold the case closed. Holds up to 10 tools, one ruler holder, and one Velcro sealed pouch. **Instruments not included.**

36 x 23 cm (unfolded)

Black

Vinyl

No. 20831-10



### Instrument Tip Protectors

These tip protectors are a convenient, color-coded way to protect surgical instrument tips and edges. They provide exceptional protection and are vented to prevent rusting from condensation. Instrument tip protectors are recommended during handling, storage and sterilization. **Autoclavable.**



1.6 x 19.1 mm  
Pkg. of 25  
Material: Plastisol

No. 29060-01



2 x 19 mm  
Pkg. of 25  
Material: Plastisol

No. 29060-02



2.8 x 19.1 mm  
Pkg. of 25  
Material: Plastisol

No. 29060-03



3.2 x 25.4 mm  
Pkg. of 25  
Material: Plastisol

No. 29060-04



29060-00




4.8 x 25.4 mm  
Pkg. of 25  
Material: Plastisol

No. 29060-05



9.5 x 25.4 mm  
Pkg. of 25  
Material: Plastisol

No. 29060-06



Assorted pack, five of each size  
Pkg. of 30  
Material: Plastisol

No. 29060-00

### Double Tip Protectors

These vented tip protectors are designed to hold hemostats and rongeurs firmly while allowing the jaws, box hinge and tips to stay open for thorough sterilization. **Autoclavable.**



12.7 x 2 mm  
Pkg. of 10  
Material: Plastisol

No. 29066-01



9 x 25.5 mm  
Pkg. of 10  
Material: Plastisol

No. 29066-02



15 x 25 mm  
Pkg. of 10  
Material: Plastisol

No. 29066-03



25 x 26 mm  
Pkg. of 10  
Material: Plastisol

No. 29066-04



Assorted pack,  
three of each size  
Pkg. of 12  
Material: Plastisol

No. 29066-00