

**Stainless Steel Burrs**

Length: 44 mm  
 Shaft Ø: 2.3 mm  
 Pkg. of 10 burrs

Tip Ø (mm)	Item
0.7 *	<b>No. 19008-07</b>
0.9 **	<b>No. 19008-09</b>
1.4	<b>No. 19008-14</b>

\* Use this size for screws No. 19010-10 (page 64).  
 \*\* Use this size for screws No. 19010-00 (page 64).



**Carbon Steel Burrs**

Length: 44 mm  
 Shaft Ø: 2.3 mm  
 Pkg. of 10 burrs



19007-07

Tip Ø (mm)	Item
0.5	<b>No. 19007-05</b>
0.7 *	<b>No. 19007-07</b>
0.9 **	<b>No. 19007-09</b>
1.4	<b>No. 19007-14</b>
2.1	<b>No. 19007-21</b>
2.9	<b>No. 19007-29</b>

\* Use this size for screws No. 19010-10 (page 64).  
 \*\* Use this size for screws No. 19010-00 (page 64).



**Tungsten Carbide Burrs**

Length: 44 mm  
 Shaft Ø: 2.3 mm  
 Pkg. of 5 burrs



Tip Ø (mm)	Item
0.5	<b>No. 19009-05</b>
0.7 *	<b>No. 19009-07</b>
0.9 **	<b>No. 19009-09</b>

\* Use this size for screws No. 19010-10 (page 64).  
 \*\* Use this size for screws No. 19010-00 (page 64).



**Stainless Steel Trephines**

Side opening for clearing of material.

Length: 44 mm  
 Shaft Ø: 2.3 mm  
 Sold individually

1.8 - 2.7 mm    5 mm

Tip Outer Ø (mm)	Item
1.8	<b>No. 18004-18</b>
2.3	<b>No. 18004-23</b>
2.7	<b>No. 18004-27</b>
5	<b>No. 18004-50</b>



**Stainless Steel Circular Saw**

For cutting bone and other hard substances. Attaches to mandrel (No. 18000-18), on the right side.



Ø: 1.6 cm  
 Thickness: 0.1 mm  
 Pkg. of 10 saws

**No. 18000-20**

**Stainless Steel Mandrel**

Screw allows attachment of discs and circular saws for use with high speed drills.



Length: 45 mm  
 Shaft Ø: 2.3 mm

**No. 18000-18**