

**High Temperature Cautery Kit**

This high temperature cauterizer is a great alternative to disposable cautery units. Temperature at tips is approximately 1000°C (1832°F). Replacement cautery tips are available. **Not autoclavable.**



Plastic body  
 Stainless steel tips: 0.2 mm Ø  
 Operation temperature: 1000°C (1832°F)  
 Total weight: 70 g  
 Length: 16 cm (without tip)

Cautery kit includes:

- 1 Cautery handle
- 2 AA batteries
- 2 Standard tips
- 1 Storage case
- 2 Fine tips

65% of actual size

**Low Temperature Cautery Kit**

This low temperature cautery kit is a great alternative to disposable cautery units. Temperature at the tips is approximately 649°C (1200°F). Replacement cautery tips are available. **Not autoclavable.**



Plastic body  
 Stainless steel tips: 0.2 mm Ø  
 Operation temperature: 649°C (1200°F)  
 Total weight: 60 g  
 Length: 12 cm (without tip)

Cautery kit includes:

- 1 Cautery handle
- 1 AA battery
- 2 Fine tips
- 1 Storage case
- 2 Elongated fine tips



Replacement tips  
**Single use**  
 Straight / Standard  
 Sterile: Yes

Box of 10 tips

No. 18010-00



Replacement tips  
**Single use**  
 Straight / Fine  
 Sterile: Yes

Box of 10 tips

No. 18010-02



Replacement tips  
**Single use**  
 Straight / Fine  
 Sterile: Yes

Box of 10 tips

No. 18019-01



Replacement tips  
**Single use**  
 Elongated Straight / Fine  
 Sterile: Yes

Box of 10 tips

No. 18019-02