

**Neuro Punches**

N.I.H. style neuro punches for punch removal of frozen brain microregions.



ID (mm)	OD (mm)	Gauge	Length (cm)	Item
0.31	0.61	24	4.2	<b>No. 18036-24</b>
0.41	0.7	22	4.5	<b>No. 18036-22</b>
0.69	1	19	4.5	<b>No. 18036-19</b>

**Absorption Spears and Triangles**

These synthetic surgical sponges absorb fluids twice as fast as the best cotton sponges and have three times their fluid-holding capacity. They are made from strong, soft, lint-free biocompatible-proven material that is nonirritating to delicate tissue and leaves no fibrous residue that might cause wound contamination. The Absorption Spears are sterilized by radiation.



Absorption Spears  
7 cm  
Sterile: Yes  
Pkg. of 200 spears

**No. 18105-01**



Absorption Triangles  
2 cm  
Sterile: Yes  
Pkg. of 200 triangles


**No. 18105-04**



**Sample Corers**

For quantitative analysis. Comes with plunger.

8.5 cm



ID (mm)	OD (mm)	Gauge	Item
0.35	0.8	21	<b>No. 18035-35</b>
0.5	1	19	<b>No. 18035-50</b>
0.8	1.45	17	<b>No. 18035-80</b>
1	1.8	15	<b>No. 18035-01</b>
2	3	11	<b>No. 18035-02</b>
3	4	8	<b>No. 18035-03</b>
5	6	4	<b>No. 18035-05</b>

**Fiber Optic Swabs**



Straight

15 cm  
Tip Ø: 2.5 mm  
Pack of 100

**No. 18106-01**

**Micro Cleaning Brush**



7 cm  
Nylon  
Bag of 25

**No. 18196-10**