

Vannas Scissors with Microserrations

Microserrated to minimize slippage.



Straight Sharp/Sharp

8.5 cm
Cutting edge: 7 mm
Tip Ø: 0.1 mm

No. 15007-08



Straight Sharp/Sharp

8.5 cm
Cutting edge: 5 mm
Tip Ø: 0.1 mm

No. 15070-08



Curved Sharp/Sharp

8.5 cm
Cutting edge: 5 mm
Tip Ø: 0.1 mm

No. 15071-08



Straight Sharp/Sharp

8 cm
Cutting edge: 5 mm
Tip Ø: 0.125 mm

No. 15044-08



Curved Sharp/Sharp

8 cm
Cutting edge: 5 mm
Tip Ø: 0.125 mm

No. 15043-08

Vannas



Vannas

Ultra long cutting edge without sacrificing tip width.



Straight Sharp/Sharp

8.5 cm
Cutting edge: 8 mm
Tip Ø: 0.1 mm

No. 15009-08