

**Long Handle Spring Scissors**



**Curved  
Sharp/Sharp**

18.5 cm  
Cutting edge: 5 mm  
Tip Ø: 0.2 mm

**No. 15027-18**

**Trident II**

Special design with short handle and inner spring.



**Straight  
Sharp/Sharp**

8 cm  
Cutting edge: 8 mm  
Tip Ø: 0.2 mm

**No. 15110-08**



**Angled to side  
Sharp/Sharp**

8 cm  
Cutting edge: 8 mm  
Tip Ø: 0.2 mm

**No. 15111-08**