

**Micro Scissors**

With microserrations to minimize slippage.




300%



**Straight  
Sharp/Sharp**

14.5 cm  
Cutting edge: 2.5 mm  
Tip Ø: 0.175 mm

**No. 15030-14**



300%




**Angled to side  
Sharp/Sharp**

14.5 cm  
Cutting edge: 2.5 mm  
Tip Ø: 0.175 mm

**No. 15031-14**

**Micro Scissors**

The finest blades on spring scissors. Designed for ophthalmic surgery (vitrectomy). The spring handles are made of ribbed plastic to ensure a solid grip. **Autoclavable.**



300%



**Straight  
Sharp/Sharp**

14 cm  
Cutting edge: 2.2 mm  
Tip Ø: 0.15 mm

**No. 15036-14**



300%



**Angled to side 90°  
Sharp/Sharp**

14 cm  
Cutting edge: 1 mm  
Tip Ø: 0.1 mm

**No. 15037-14**