

**Bonn-Gawenda Scissors**

The innovative Bonn-Gawenda scissors from our Exclusive Line combine a precise cutting edge, similar to our spring scissors and our ToughCut® feature which ensures an exact cut without slipping. These scissors have perfect tips for precision work under a microscope. Designed with comfort and ergonomics in mind, the large finger loops and specially engineered handle ensure exceptional control throughout each use.










  
**Straight Sharp/Sharp**
  
 10.5 cm
   
 Cutting edge: 10 mm
   
**No. 14854-10**



  
**Curved Sharp/Sharp**
  
 10.5 cm
   
 Cutting edge: 10 mm
   
**No. 14855-10**


**Extra Fine Bonn**


Our finest iris scissors, best suited for microdissection under high magnification.






  
**Straight Sharp/Sharp**
  
 8.5 cm
   
 Cutting edge: 10 mm
   
**No. 14084-08**


  
**Curved Sharp/Sharp**
  
 8.5 cm
   
 Cutting edge: 10 mm
   
**No. 14085-08**


  
**Straight Sharp/Blunt**
  
 8.5 cm
   
 Cutting edge: 10 mm
   
**No. 14083-08**