

**Extra Fine Bonn**

Highly polished.



**ToughCut® Bonn-Strabismus**

A fusion of two light and delicate iris scissors with ToughCut cutting edges.



**Hardened Fine Iris Scissors**

These scissors have an anti-reflective finish to reduce glare.



**Straight Sharp/Sharp**

8.5 cm  
Cutting edge: 10 mm

**No. 14184-09**



**Curved Sharp/Sharp**

8.5 cm  
Cutting edge: 10 mm

**No. 14185-09**



**Straight Blunt/Blunt**

8.5 cm  
Cutting edge: 10 mm

**No. 14084-09**



**Curved Blunt/Blunt**

8.5 cm  
Cutting edge: 10 mm

**No. 14085-09**



**Straight Sharp/Sharp**

9 cm  
Cutting edge: 16 mm

**No. 14090-09**

11 cm  
Cutting edge: 16 mm

**No. 14090-11**



**Curved Sharp/Sharp**

9 cm  
Cutting edge: 16 mm

**No. 14091-09**

11 cm  
Cutting edge: 16 mm

**No. 14091-11**