


Dissector Scissors

Slim blades for cutting tissue.

**Angled to side
Sharp/Sharp**10 cm
Cutting edge: 6 mm**No. 14081-09****Dissector Scissors**

Similar to 14081-09, but with stronger blades for cutting fine bone.

**Angled to side
Sharp/Sharp**10.5 cm
Cutting edge: 8 mm**No. 14082-09****Spencer Ligature**

Used for cutting stitches or suture.

**Straight
Sharp/Sharp**9 cm
Cutting edge: 16 mm**No. 14076-09**