

Micro Friedmann-Pearson Rongeur

Experience the unparalleled precision of our flagship rongeur which features a specialized 0.5 mm cup diameter, tailored for the finest research techniques. These rongeurs guarantee effortless cutting while ensuring captured tissue remains secure and intact.



200%




CAUTIONARY NOTE
Cautionary note: The fine rongeur tips shown on this page are extremely fragile. Removal of material with a twisting motion can damage the tips.



Straight

14 cm
Cup size: 0.5 mm

No. 16220-14



Curved

14 cm
Cup size: 0.5 mm

No. 16221-14