



Friedmann-Pearson Rongeur

Best suited for delicate work, this is one of our finest tip rongeurs in our product line. The fine-cupped jaws have been designed with a longer and slimmer approach for better visibility.




Cautionary note: The fine rongeur tips shown on this page are extremely fragile. Removal of material with a twisting motion can damage the tips.



Straight

14 cm
Cup size: 0.7 mm

No. 16120-14



Curved

14 cm
Cup size: 0.7 mm


No. 16121-14



Straight

14 cm
Cup size: 1 mm

No. 16020-14



Curved

14 cm
Cup size: 1 mm

No. 16021-14