


Micro Curette

Curette tips consist of a scoop with sharp edges and are used for scraping or cleaning. They can also be used for separating and transporting, inoculating plates or transferring small amounts of chemicals.




Serrations on cutting edge




13 cm

No. 10080-05




13 cm

No. 10081-10



13 cm

No. 10082-15




13 cm

No. 10083-20

Delicate Bone Scraper


For scraping tissue without cutting into bone.



16.5 cm

No. 10075-16


Austin Chisel



16 cm

No. 10096-16

Dissecting Chisel



12.5 cm

No. 10095-12