

Scalpel Handle #4 with Silicone

For heavy, gross anatomy work.
Molded to fit hand.



13.5 cm

No. 10009-13

Standard Scalpel Handle #4 Solid

Especially suitable for large dissections.



13 cm

No. 10004-13

Blades for #4 Scalpel Handle

The alloy used has a high carbon content which provides a longer lasting edge. Blades come individually wrapped. **Non-sterile.**



#20 blades

Box of 100 blades

No. 10020-00



#21 blades

Box of 100 blades

No. 10021-00



#22 blades

Box of 100 blades

No. 10022-00



#23 blades

Box of 100 blades

No. 10023-00